

**NAME**

d2scan – two-motor scan relative to the starting positions

**SYNOPSIS**

```
d2scan motor1 start1 end1 motor2 start2 end2 intervals time
```

**DESCRIPTION**

d2scan scans two motors, as specified by *motor1* and *motor2*. Each motor moves the same number of intervals. If each motor is at a position  $X$  before the scan begins, it will be scanned from  $X+start$  to  $X+end$ . The step size for each motor is  $(start - end) / intervals$ . The number of data points collected will be  $intervals + 1$ . Count time is given by *time*, which if positive, specifies seconds and if negative, specifies monitor counts.

Upon termination, the motors are returned to their starting positions.

**EXAMPLE**

```
d2scan th -1 1 phi -2 2 20 10
```

**SEE ALSO**

motors ascan a2scan a3scan dscan d3scan mesh